

# The Relaxatio Sales Pitch

I would like to present to you the newest technology currently available for pain management, which will give you immediate results.

This is Relaxatio. It uses a powerful combination of red and infrared LED lights, which provide a deeper and more effective treatment for all skin types. Relaxatio was designed to eliminate pain.

Pain is normally caused by damage to our tissues. As our body grows older, it stores less energy. This decrease in energy results in longer recovery time from injuries and higher levels of inflammation, which can lead to chronic pain.

The Relaxatio LED light increases ATP, which is stored energy in our body. This increased energy accelerates the repair processes of the cells and reduces inflammation. Additionally, the process stimulates blood circulation to the treated area, bringing more proteins and oxygen that are essential for recovery and healing.

The red and infrared lights penetrate 3mm (millimeter) under the skin, treating all layers of the skin from surface to core — epidermis, dermis, all the way down to the innermost and thickest layer of the skin, the hypodermis. Relaxatio comes with a glass treatment probe for optimal light exposure. The glass probe heats-up to 42° Celsius (107° Fahrenheit) to provide topical heating, which further increases local blood flow at the treated area, providing extra nutrients and extra oxygen for better healing of the skin.

The FDA-approved Relaxatio device is clinically proven and was found to be completely safe to use and with no reported side effects. In order to treat pain, it is recommended to first use the Relaxatio device by itself on clean skin, for up to 10 minutes per treated area, and then add the Red Algae Infused X Pain Relief Gel (provided) for five more minutes of treatment with the device. Treat each painful area as required to eliminate the pain. For increased effectiveness, use the X Pain Relief Gel after treatment as well. Red algae is a well-known miracle ingredient, with incredible therapeutic benefits for the skin. It conditions, hydrates, softens and detoxifies the skin, replenishing essential vitamins and minerals, as well as reducing pain and the effects of inflammation. LED light therapy is used today by the world's leading dermatologists due to its proven benefits for the skin. It is currently being used by medical doctors all over the world who are treating chronic joint pain, inflammation, nerve damage, and athletic recovery.



U.S. Patent No.  
D 16/013,068



# Treatment - simple, fast and safe

1. Make sure the Relaxatio is fully charged before use (green light LED indicator).
2. After you turn on the device, it will start heating up slowly (no need to wait). Place the glass probe on your skin and the red light will automatically activate.
3. Apply light pressure while slowly moving the Relaxatio in small circular motions.
4. Treat the painful areas: lower back, neck, shoulder, etc.
5. It is recommended to use the Relaxatio with a small amount of pain relief cream or gel (Red Algae Infused X Pain Relief Gel) on the treated areas during each session and to massage it gently in circular motions. Any cream you use with the treatment will make it more effective.
6. Clean the Relaxatio after every use with a cloth, cotton ball or paper towel.

## Relaxatio™

